

# Penetration Test Report

## CardioSense SaMD Application Security Assessment

---

Report Reference:	PT-2026-001
Test Date:	2026-01-08 to 2026-01-12
Testing Firm:	SecureHealth Consulting Ltd
Lead Tester:	Dr. Marcus Webb, OSCP, CISSP
Application:	CardioSense v2.1.0 (SaMD Class C)
Scope:	Web API, Mobile App, Cloud Infrastructure
Classification:	Confidential - Restricted Distribution

### 1. Executive Summary

SecureHealth Consulting performed a comprehensive penetration test of the CardioSense SaMD application. Testing covered OWASP Top 10, API security, authentication/authorization, data encryption, and cloud infrastructure configuration. Testing was performed against the staging environment which mirrors production configuration.

Overall Risk Rating: LOW. No critical or high-severity vulnerabilities were identified. Three medium-severity and five low-severity findings were documented.

### 2. Findings Summary

Severity	ID	Description
MEDIUM	PT-F001	Missing rate limiting on /api/auth/login endpoint
MEDIUM	PT-F002	Session token not rotated after privilege change
MEDIUM	PT-F003	Verbose error messages expose stack traces in staging
LOW	PT-F004	Missing Content-Security-Policy header on web app
LOW	PT-F005	X-Powered-By header reveals server technology
LOW	PT-F006	Autocomplete not disabled on sensitive form fields
LOW	PT-F007	Cookie missing SameSite=Strict attribute
LOW	PT-F008	DNS CAA records not configured

### 3. Detailed Findings

#### **PT-F001: Missing Rate Limiting (Medium)**

The /api/auth/login endpoint does not enforce rate limiting, allowing an attacker to perform brute-force attacks. Recommendation: Implement rate limiting of 5 attempts per minute with exponential backoff.

#### **PT-F002: Session Token Rotation (Medium)**

After a user changes their role or permissions, the session token is not regenerated. Recommendation: Rotate session tokens after any privilege change to prevent session fixation attacks.

#### **PT-F003: Verbose Error Messages (Medium)**

API error responses include full stack traces in the staging environment. Recommendation: Ensure production builds strip stack traces from error responses. Verify this is not the case in production deployment.

### 4. Conclusion & Recommendations

The CardioSense application demonstrates a strong security posture overall. No critical vulnerabilities were identified that would allow unauthorized access to patient data or medical device functionality. The medium-severity findings should be addressed within 30 days. Low-severity findings should be addressed within the next release cycle.

The application meets the security requirements specified in IEC 62443 and the FDA Guidance on Cybersecurity for Medical Devices (2023).

---

**Dr. Marcus Webb, OSCP, CISSP**

Principal Security Consultant

SecureHealth Consulting Ltd

Date: 2026-01-15